

## Product Datasheet

### Benzoylecgonine (BE)-BTG Antigen LA466

<b>Product Name</b>	Benzoylecgonine (BE)-BTG Antigen
<b>Catalog Number</b>	LA466
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Antigen-BTG conjugate
<b>Storage</b>	Short Term ( $\leq 1$ week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	$\geq 90\%$ by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Matched Pairs</b>	HM777 (detection)
<b>Preservatives</b>	0.01% Sodium Azide