

Product Datasheet

Anti-Mullerian Hormone (AMH) Antigen, Recombinant >90% LA467

Product Name	Anti-Mullerian Hormone (AMH) Antigen, Recombinant >90%
Catalog Number	LA467
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Recombinant Antigen
Storage	Short Term (\leq 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	>90%
Buffer	20 mM Phosphate with 0.01% SDS
Molecular Weight	45 kDa
Matched Pairs	HM733 or HM737
Expression System	E. coli