

## Product Datasheet

### Astrovirus Antigen, Recombinant LA470

**Product Name** Astrovirus Antigen, Recombinant

**Catalog Number** LA470

**Tested Applications** LF, ELISA

**Brand** EastCoast Bio

**Form/Appearance** Recombinant Antigen

**Storage** Short Term ( $\leq 2$  weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Buffer** TBS, pH 8.0

**Matched Pairs** HM803, HM804

**Expression System** E. coli

**Preservatives** 0.09% Sodium Azide