

Product Datasheet

Astrovirus Antigen, Recombinant LA470

Product Name Astrovirus Antigen, Recombinant

Catalog Number LA470

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Recombinant Antigen

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Buffer TBS, pH 8.0

Matched Pairs HM803, HM804

Expression System E. coli

Preservatives 0.09% Sodium Azide