

Product Datasheet

S100A12 Antigen, Recombinant LA479

Product Name	S100A12 Antigen, Recombinant
Catalog Number	LA479
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Recombinant Antigen
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	TBS, pH 7.4
Matched Pairs	JP135, JP136
Expression System	E. coli
Preservatives	0.09% Sodium Azide