

Product Datasheet

Elastase CELA3A Antigen, Recombinant LA480

Product Name	Elastase CELA3A Antigen, Recombinant
Catalog Number	LA480
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Recombinant Antigen
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	≥ 95 %
Buffer	Tris buffer with additives and salts, pH 8.0
Matched Pairs	JP137, JP138
Expression System	Mammalian Cell
Preservatives	0.09% Sodium Azide