

## Product Datasheet

### D-dimer Antigen, Native ≥90% LA482

<b>Product Name</b>	D-dimer Antigen, Native ≥90%
<b>Catalog Number</b>	LA482
<b>Tested Applications</b>	Control
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Native Antigen
<b>Storage</b>	Short Term (≤ 1 week): 2-8°C. Long Term: -70°C. Avoid repeated freezing and thawing.
<b>Purity</b>	≥90% by SDS-PAGE
<b>Source</b>	Human Plasma
<b>Buffer</b>	20 mM Tris-HCl, 50 mM NaCl, 2 mM KCl, pH 7.4
<b>Molecular Weight</b>	180 kDa
<b>Matched Pairs</b>	HM812
<b>Preservatives</b>	0.1% Proclin 300