

Product Datasheet

Vitamin B 12-BSA Antigen LA484

Product Name	Vitamin B 12-BSA Antigen
Catalog Number	LA484
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM830, HM831 or HM832 (detection)
Preservatives	0.01% Sodium Azide