

Product Datasheet

C-Reactive Protein, Canine (cCRP) Antigen, Native LA485

Product Name	C-Reactive Protein, Canine (cCRP) Antigen, Native
Catalog Number	LA485
Tested Applications	LF, ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Lyophilized; Purified Native Antigen
Storage	Lyophilized: 2-8°C. Reconstituted: 2-8°C < 1 week; long term storage @ -20°C or below. Avoid repeated freezing and thawing.
Note	Reconstitute with 1 ml of water. Aliquot to use size prior to freezing.
Source	Dog serum
Preservatives	0.05% Proclin 300