

Product Datasheet

Ebola Virus (Sudan-NP) Antigen, Recombinant >95% LA487

Product Name Ebola Virus (Sudan-NP) Antigen, Recombinant >95%

Catalog Number LA487

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen

Storage -20°C. Avoid repeated freezing and thawing.

Purity >95% by SDS-PAGE

Buffer 10 mM Phosphate Buffered Saline, pH 7.2

Specificity Ebola Virus Sudan nucleoprotein

Molecular Weight 15 kDa

Expression System E.coli

Preservatives 0.02% Sodium Azide