

## Product Datasheet

### **Ebola Virus (Sudan-NP) Antigen, Recombinant >95% LA487**

**Product Name** Ebola Virus (Sudan-NP) Antigen, Recombinant >95%

**Catalog Number** LA487

**Tested Applications** LF, ELISA

**Brand** EastCoast Bio

**Form/Appearance** Purified Recombinant Antigen

**Storage** -20°C. Avoid repeated freezing and thawing.

**Purity** >95% by SDS-PAGE

**Buffer** 10 mM Phosphate Buffered Saline, pH 7.2

**Specificity** Ebola Virus Sudan nucleoprotein

**Molecular Weight** 15 kDa

**Expression System** E.coli

**Preservatives** 0.02% Sodium Azide