

## Product Datasheet

### **Ebola Virus (Zaire-gp) Antigen, Recombinant >95% LA488**

<b>Product Name</b>	Ebola Virus (Zaire-gp) Antigen, Recombinant >95%
<b>Catalog Number</b>	LA488
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen (His-tag)
<b>Storage</b>	-20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	>95% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2 with 25 mM Arginin
<b>Specificity</b>	Ebola Virus Zaire glycoprotein
<b>Molecular Weight</b>	38 kDa
<b>Expression System</b>	E.coli
<b>Preservatives</b>	0.02% Sodium Azide