

Product Datasheet

Ebola Virus (Zaire-gp) Antigen, Recombinant >95% LA488

Product Name	Ebola Virus (Zaire-gp) Antigen, Recombinant >95%
Catalog Number	LA488
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 with 25 mM Arginin
Specificity	Ebola Virus Zaire glycoprotein
Molecular Weight	38 kDa
Expression System	E.coli
Preservatives	0.02% Sodium Azide