

## Product Datasheet

### Inhibin B (INHB) Antigen, Recombinant LA496

<b>Product Name</b>	Inhibin B (INHB) Antigen, Recombinant
<b>Catalog Number</b>	LA496
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen (His-tag)
<b>Storage</b>	Short Term: 2-8°C. Long Term: ≤-70°C
<b>Purity</b>	>90%
<b>Buffer</b>	20 mM Phosphate Buffer, pH 7.4 with 0.01% SDS
<b>Matched Pairs</b>	HM909, HM910
<b>Expression System</b>	E. coli