

Product Datasheet

Inhibin B (INHB) Antigen, Recombinant LA496

Product Name	Inhibin B (INHB) Antigen, Recombinant
Catalog Number	LA496
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	Short Term: 2-8°C. Long Term: ≤-70°C
Purity	>90%
Buffer	20 mM Phosphate Buffer, pH 7.4 with 0.01% SDS
Matched Pairs	HM909, HM910
Expression System	E. coli