

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1-gp41) Antigen, Recombinant >90% LA497

Product Name	Human Immunodeficiency Virus 1 (HIV 1-gp41) Antigen, Recombinant >90%
Catalog Number	LA497
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	2-8°C. Do not freeze.
Purity	>90%
Buffer	20 mM Hepes, 300 mM NaCl, 102 mM Imidazole, 0.2% SDS
Specificity	gp41
Molecular Weight	18 kDa
Expression System	E. coli
Preservatives	0.05% Proclin 300