

## Product Datasheet

### Human Immunodeficiency Virus 2 (HIV 2-gp36) Antigen, Recombinant >90% LA498

<b>Product Name</b>	Human Immunodeficiency Virus 2 (HIV 2-gp36) Antigen, Recombinant >90%
<b>Catalog Number</b>	LA498
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen (His-tag)
<b>Storage</b>	2-8°C. Do not freeze.
<b>Purity</b>	>90%
<b>Buffer</b>	20 mM HEPES, 300 mM NaCl, 0.2% SDS
<b>Specificity</b>	gp36
<b>Molecular Weight</b>	20 kDa
<b>Expression System</b>	E. coli
<b>Preservatives</b>	0.05% Proclin 300