

Product Datasheet

Hepatitis B Surface Antigen (HBsAg-subtype Ay), Native >99% LA501

Product Name Hepatitis B Surface Antigen (HBsAg-subtype Ay), Native >99%

Catalog Number LA501

Tested Applications ELISA

Brand EastCoast Bio

Form/Appearance Purified & 0.2 µm filtered Native Antigen

Storage 2-8°C

Note Tested negative for HIV 1/2 and HCV, however as with all biological material it should be handled as potentially infectious.

Purity >99%

Source Human Plasma

Buffer 25 mM Phosphate Buffer, pH 7.4 with 150 mM NaCl, 5% sucrose

Preservatives 0.1% Sodium Azide