
Product Datasheet

Human Chorionic Gonadotropin (HCG-intact) Antigen, Native >95% LA507L

Product Name	Human Chorionic Gonadotropin (HCG-intact) Antigen, Native >95%
Catalog Number	LA507L
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Lyophilized Native Antigen
Storage	Lyophilized: Short Term 2-8°C Long Term: -20°C. Reconstituted: -20°C. Avoid repeated freezing and thawing
Purity	>95% by SDS-PAGE
Source	Human Urine
Molecular Weight	36 kDa