

Product Datasheet

Human Chorionic Gonadotropin ? (HCG?) Antigen, Native >95% LA510L

Product Name Human Chorionic Gonadotropin β (HCG β) Antigen, Native >95%

Catalog Number LA510L

Tested Applications Control

Brand EastCoast Bio

Form/Appearance Lyophilized Native Antigen

Storage Lyophilized: Short Term 2-8°C Long Term: -20°C. Reconstituted: -20°C. Avoid repeated freezing and thawing

Purity >95% by SDS-PAGE

Source Human Urine

Molecular Weight 22 kDa