

## Product Datasheet

### Adenovirus Hexon Antigen, Recombinant ≥90% LA514

<b>Product Name</b>	Adenovirus Hexon Antigen, Recombinant ≥90%
<b>Catalog Number</b>	LA514
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	≥90%
<b>Buffer</b>	Tris buffer with additives and salts, pH 8.0
<b>Matched Pairs</b>	HM948
<b>Expression System</b>	Mammalian Cell
<b>Preservatives</b>	0.09% Sodium Azide