

Product Datasheet

Adenovirus Antigen, Native LA515

Product Name	Adenovirus Antigen, Native
Catalog Number	LA515
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Native Antigen
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Source	Mammalian Cells
Buffer	Phosphate Buffered Saline with BSA
Matched Pairs	HM948
Preservatives	0.09% Sodium Azide