

Product Datasheet

Legionella pneumophila Antigen, Native LA516

Product Name	Legionella pneumophila Antigen, Native
Catalog Number	LA516
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Inactivated Native Antigen
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Source	L. pneumophila culture
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM949, HM950
Preservatives	0.09% Sodium Azide