

Product Datasheet

S100B Antigen, Recombinant >95% LA521

Product Name	S100B Antigen, Recombinant >95%
Catalog Number	LA521
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	>95%
Buffer	Tris Buffer, pH 8.0
Matched Pairs	JP152, JP153
Expression System	E. coli
Preservatives	0.09% Sodium Azide