

## Product Datasheet

### Human Papilloma Virus (HPV-18) Antigen, Recombinant >90% LA524

**Product Name** Human Papilloma Virus (HPV-18) Antigen, Recombinant >90%

**Catalog Number** LA524

**Tested Applications** Control

**Brand** EastCoast Bio

**Form/Appearance** Purified Recombinant Antigen (GST-tag)

**Storage** Short Term ( $\leq 1$  week) 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Purity** >90% by SDS-PAGE

**Buffer** 10 mM Phosphate Buffered Saline and 2M Urea

**Molecular Weight** 78 kDa

**Expression System** E. coli