

## Product Datasheet

### **Clostridium difficile toxin A Antigen, Native >95% LA528L**

<b>Product Name</b>	Clostridium difficile toxin A Antigen, Native >95%
<b>Catalog Number</b>	LA528L
<b>Tested Applications</b>	Control
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Lyophilized Native Antigen
<b>Storage</b>	Lyophilized: 2-8°C. Reconstituted: Short Term ( $\leq$ 1 week) 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Reconstitute with distilled water. Do not vortex.
<b>Purity</b>	>95% by SDS-PAGE
<b>Source</b>	C. difficile culture