

Product Datasheet

Clostridium difficile toxoid B Antigen, Native >95% LA531L

Product Name	Clostridium difficile toxoid B Antigen, Native >95%
Catalog Number	LA531L
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Lyophilized Native Antigen
Storage	Lyophilized: 2-8°C. Reconstituted: Short Term (\leq 1 week) 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Reconstitute with distilled water. Do not vortex.
Purity	>95% by SDS-PAGE
Source	C. difficile toxoid from culture