

Product Datasheet

Hepatitis B Surface Antigen (HBsAg-subtype Ad), Recombinant >90% LA549

Product Name	Hepatitis B Surface Antigen (HBsAg-subtype Ad), Recombinant >90%
Catalog Number	LA549
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	≥90%
Source	CHO cells
Buffer	10 mM Phosphate Buffer, pH 7.4 with 150 mM NaCl
Expression System	Mammalian Cell (CHO)
Preservatives	0.1% Proclin 300