

## Product Datasheet

### Coronavirus (CoV-NP 229E) Antigen, Recombinant >95% LA557

<b>Product Name</b>	Coronavirus (CoV-NP 229E) Antigen, Recombinant >95%
<b>Catalog Number</b>	LA557
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen (His-tag)
<b>Storage</b>	Short Term ( $\leq$ 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Peptide composed of 359 amino acids containing immunodominant regions.
<b>Purity</b>	>95%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2 with 3M Urea
<b>Expression System</b>	E. coli