

## Product Datasheet

### Coronavirus (CoV-NP-NL63; 221-340) Antigen, Recombinant >95% LA558

<b>Product Name</b>	Coronavirus (CoV-NP-NL63; 221-340) Antigen, Recombinant >95%
<b>Catalog Number</b>	LA558
<b>Tested Applications</b>	ELISA, WB
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen
<b>Storage</b>	Short Term ( $\leq$ 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Peptide composed of 130 total amino acids including 221-340.
<b>Purity</b>	>95%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Expression System</b>	E. coli