

## Product Datasheet

### Coronavirus (SARS-CoV ENV; 1-76) Antigen, Recombinant >95% LA561

**Product Name** Coronavirus (SARS-CoV ENV; 1-76) Antigen, Recombinant >95%

**Catalog Number** LA561

**Tested Applications** ELISA, WB

**Brand** EastCoast Bio

**Form/Appearance** Purified Recombinant Antigen (GST-tag)

**Storage** Short Term ( $\leq 1$  week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Note** Peptide composed of amino acids 1-76 of envelope protein. MW= 34 kDa

**Purity** >95%

**Buffer** 50 mM Tris-HCl, 60 mM NaCl with 50% glycerol

**Expression System** E. coli