

Product Datasheet

Coronavirus (SARS-CoV Matrix; 182-216) Antigen, Recombinant >95% LA562

Product Name	Coronavirus (SARS-CoV Matrix; 182-216) Antigen, Recombinant >95%
Catalog Number	LA562
Tested Applications	ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (GST-tag)
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of amino acids 182-216 of matrix protein.
Purity	>95%
Buffer	50 mM Tris-HCl, 60 mM NaCl with 50% glycerol
Molecular Weight	30 kDa
Expression System	E. coli