

## Product Datasheet

### Coronavirus (SARS-CoV NP N-Terminus; 1-49, 192-220) Antigen, Recombinant >95% LA565

<b>Product Name</b>	Coronavirus (SARS-CoV NP N-Terminus; 1-49, 192-220) Antigen, Recombinant >95%
<b>Catalog Number</b>	LA565
<b>Tested Applications</b>	ELISA, WB
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen (GST-tag)
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Peptide composed of amino acids 1-49, 192-220 of nucleoprotein.
<b>Purity</b>	>95%
<b>Buffer</b>	50 mM Tris-HCl, 60 mM NaCl with 50% glycerol
<b>Molecular Weight</b>	35 kDa
<b>Expression System</b>	E. coli