

Product Datasheet

Coronavirus (SARS-CoV NP N-Terminus; 1-49) Antigen, Recombinant >95% LA566

Product Name	Coronavirus (SARS-CoV NP N-Terminus; 1-49) Antigen, Recombinant >95%
Catalog Number	LA566
Tested Applications	ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (GST-tag)
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of amino acids 1-49 of nucleoprotein.
Purity	>95%
Buffer	50 mM Tris-HCl, 60 mM NaCl with 50% glycerol
Molecular Weight	32 kDa
Expression System	E. coli