

## Product Datasheet

### Coronavirus (SARS-CoV Spike Proten S1; 12-53, 90-115, 171-203) Antigen, Recombinant >95% LA567

**Product Name** Coronavirus (SARS-CoV Spike Proten S1; 12-53, 90-115, 171-203)  
Antigen, Recombinant >95%

**Catalog Number** LA567

**Tested Applications** ELISA, WB

**Brand** EastCoast Bio

**Form/Appearance** Purified Recombinant Antigen (GST-tag)

**Storage** Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Note** Peptide composed of amino acids 12-53, 90-115, 171-203 of spike protein S1 domain.

**Purity** >95%

**Buffer** 25 mM Tris-HCl with detergents and 50% glycerol

**Molecular Weight** 38 kDa

**Expression System** E. coli