

Product Datasheet

Coronavirus (SARS-CoV Spike Protein S1; 408-470, 540-573) Antigen, Recombinant >95% LA568

Product Name	Coronavirus (SARS-CoV Spike Protein S1; 408-470, 540-573) Antigen, Recombinant >95%
Catalog Number	LA568
Tested Applications	ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (GST-tag)
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of amino acids 408-470, 540-573 of spike protein S2 domain.
Purity	>95%
Buffer	25 mM Tris-HCl with detergents and 50% glycerol
Molecular Weight	38 kDa
Expression System	E. coli