

Product Datasheet

Coronavirus (SARS-CoV ENV; 1-76) Antigen, Recombinant >90% LA570

Product Name Coronavirus (SARS-CoV ENV; 1-76) Antigen, Recombinant >90%

Catalog Number LA570

Tested Applications ELISA

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen (His-tag)

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Note Peptide composed of amino acids 1-76 of envelope protein.

Purity >90%

Buffer 10 mM Phosphate Buffered Saline, pH 7.2

Expression System E.coli