

## Product Datasheet

### Coronavirus (SARS-CoV Spike Protein S1; 1-53, 90-115, 171-205) Antigen, Recombinant >90% LA574

<b>Product Name</b>	Coronavirus (SARS-CoV Spike Protein S1; 1-53, 90-115, 171-205) Antigen, Recombinant >90%
<b>Catalog Number</b>	LA574
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen (His-tag)
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Peptide composed of amino acids 1-53, 90-115, 171-205 of spike protein S1 domain.
<b>Purity</b>	>90%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Expression System</b>	E.coli