

Product Datasheet

Coronavirus (SARS-CoV Spike Protein S1; 408-470, 540-573) Antigen, Recombinant >90% LA575

Product Name	Coronavirus (SARS-CoV Spike Protein S1; 408-470, 540-573) Antigen, Recombinant >90%
Catalog Number	LA575
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of amino acids 408-470, 540-573 of spike protein S2 domain.
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Expression System	E.coli