

Product Datasheet

Coronavirus (SARS-CoV Spike Protein; 20-210) Antigen, Recombinant >90% LA576

Product Name	Coronavirus (SARS-CoV Spike Protein; 20-210) Antigen, Recombinant >90%
Catalog Number	LA576
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of amino acids 20-210 of spike protein.
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Molecular Weight	32KDa
Expression System	E.coli