

Product Datasheet

Hepatitis Be Antigen (HBeAg VLP), Recombinant >95% LA577

Product Name	Hepatitis Be Antigen (HBeAg VLP), Recombinant >95%
Catalog Number	LA577
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Recombinant Virus Like Particle Antigen, (His-tag)
Storage	Short Term (≤ 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Molecular Weight	18 kDa
Expression System	E. coli
Preservatives	0.05% Sodium Azide