

## Product Datasheet

### Coronavirus (COVID-19 Spike Protein S1; 1-674) Antigen, Recombinant >95% LA582

<b>Product Name</b>	Coronavirus (COVID-19 Spike Protein S1; 1-674) Antigen, Recombinant >95%
<b>Catalog Number</b>	LA582
<b>Tested Applications</b>	Control
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein G Purified Recombinant Antigen (sFc-tag)
<b>Storage</b>	-70°C. Avoid repeated freezing and thawing
<b>Note</b>	Peptide composed of amino acids 1-674 of spike protein S1 domain.
<b>Purity</b>	>95%
<b>Buffer</b>	DPBS, pH 7.4
<b>Specificity</b>	Wuhan-Hu-1
<b>Expression System</b>	Mammalian Cell (HEK293)