

Product Datasheet

Coronavirus (COVID-19 Spike Protein S2; 685-1211) Antigen, Recombinant >80% LA583

Product Name	Coronavirus (COVID-19 Spike Protein S2; 685-1211) Antigen, Recombinant >80%
Catalog Number	LA583
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Recombinant Antigen (sFc-tag)
Storage	-70°C. Avoid repeated freezing and thawing
Note	Peptide composed of amino acids 685-1211 of spike protein S2 domain.
Purity	>80%
Buffer	DPBS, pH 7.4
Specificity	Wuhan-Hu-1
Expression System	Mammalian Cell (HEK293)