

Product Datasheet

Human IL-12p40/p70 ELISA Set 851.880

Product Name	Human IL-12p40/p70 ELISA Set
Catalog Number	851.880
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	62.5-2000 pg/ml
Sensitivity	20 pg/ml
Tested Applications	ELISA
Alternative Names	IL-12/IL23p40
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes the p40 monomer, dimer or heterodimer (p40+p35) forms of the molecule
Sample Types	Serum, Plasma, Cell culture supernatant, buffered solution