
Product Datasheet

Methylone-BSA Antigen LA608

Product Name	Methylone-BSA Antigen
Catalog Number	LA608
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM1067 (detection)
Preservatives	0.01% Sodium Azide