

Product Datasheet

Coronavirus (COVID-19 Spike Protein RBD; 331-524) Antigen, Recombinant >95% LA609

Product Name	Coronavirus (COVID-19 Spike Protein RBD; 331-524) Antigen, Recombinant >95%
Catalog Number	LA609
Tested Applications	LF, ELISA, Control
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (mFc-tag)
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of amino acids 331-524 of spike protein receptor binding domain (RBD).
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Molecular Weight	51 kDa
Expression System	Mammalian Cell (HEK293)