

## Product Datasheet

### Coronavirus (COVID-19 Spike Protein S1; 1-681) Antigen, Recombinant >95% LA610

<b>Product Name</b>	Coronavirus (COVID-19 Spike Protein S1; 1-681) Antigen, Recombinant >95%
<b>Catalog Number</b>	LA610
<b>Tested Applications</b>	LF, ELISA, Control
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen (His-tag)
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Peptide composed of amino acids 1-681 of spike protein S1 domain.
<b>Purity</b>	>95%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.0
<b>Molecular Weight</b>	76 kDa
<b>Expression System</b>	Mammalian Cell (HEK293)
<b>Preservatives</b>	0.09% Sodium Azide