

Product Datasheet

Coronavirus (COVID-19 NP; 1-419) Antigen, Recombinant >90% LA612

Product Name Coronavirus (COVID-19 NP; 1-419) Antigen, Recombinant >90%

Catalog Number LA612

Tested Applications LF

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen (His-tag)

Storage Short Term (<2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing

Note Peptide composed of amino acids 1-419 of nucleoprotein.

Purity >90%

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Molecular Weight 48 kDa

Expression System E. coli