

Product Datasheet

Coronavirus (COVID-19 NP; 1-419) Antigen, Recombinant >90% LA616

Product Name Coronavirus (COVID-19 NP; 1-419) Antigen, Recombinant >90%

Catalog Number LA616

Tested Applications ELISA

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen

Storage Short Term (<2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing

Note Peptide composed of amino acids 1-419 of nucleoprotein.

Purity >90%

Buffer 25 mM Tris-HCl, pH 8.0

Molecular Weight 72 kDa

Expression System E. coli