

## Product Datasheet

### Treponema pallidum (TpN15, TpN17, TpN47) Antigen, Recombinant >90% LA618

<b>Product Name</b>	Treponema pallidum (TpN15, TpN17, TpN47) Antigen, Recombinant >90%
<b>Catalog Number</b>	LA618
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	Syphilis
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen
<b>Storage</b>	2-8°C
<b>Purity</b>	>90%
<b>Source</b>	
<b>Buffer</b>	20 mM Tris-HCl, pH 8.0 with 5 mM EDTA and 0.1% SDS
<b>Specificity</b>	TpN15, TpN17, TpN47
<b>Matched Pairs</b>	LA619 (detection)
<b>Expression System</b>	E. coli